


6000rpm Low Speed Blood Centrifuge Machine For Laboratory With Airtight Rotors

Basic Information

- Place of Origin: CHINA
- Brand Name: MKE
- Certification: CE, ISO13485,ISO9001, SFDA
- Model Number: MK540
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Wooden Box
- Delivery Time: 7-15 work days
- Payment Terms: L/C,D/A,D/P,T/T,Western Union,MoneyGram
- Supply Ability: 100 set/month





Product Specification

- Applications: Industrial Labs,Hospitals, Veterinary Clinics, Food Labs, Universities & Colleges, Pharmaceutical Labs, Cell Culture Labs
- Max.Speed: 6000rpm
- Max.RCF: 6080xg
- Max.Capacity: 4*750ml
- Matched Rotors: Fixed Angle Rotor, Square Swing Bucket, Round Swing Bucket With Air Tight Cover
- Display: LCD
- Warranty: 1 Year Warranty For The Whole Machine, 3 Years Warranty For The Compressor And 5 Years Warranty For The Motor. Free Replacement Parts And Shipping Within Warranty.
- Highlight: **Laboratory blood bank centrifuge machine, LCD blood bank centrifuge machine, LCD blood bank centrifuge**



More Images



Product Description

The large centrifuge MK540 is of bio-safety type with good aerosols filtering system, built-in ultraviolet sterilization lamp and separate air tight cover for swing bucket.

Design

- LCD display. The centrifuge parameters can be set in operation. User-friendly program design. The adjustable modes of short spin, timing and rotational speed make the test more efficient and stable.
- Impulse centrifugation. Instant speed up to the setting value. Auto brake and fast shutdown.
- DC brushless motor with free maintenance. Faster speed acceleration and deceleration. More samples centrifugation in short time.



User friendly


- Unique appearance design with small size. Easy to use. The lid can open 90° to handle the microplates.
- Disperse the hanging liquid on tube wall quickly before or after PCR test, which improves test effects greatly.



Safety

- Produced according to international safety regulations (e.g. IEC 61010)
- Safety Electrical lid interlock, door can not open when rotor is running and machine won't run when lid open and error code appears.

- Emergency lid-lock release when power failure.
- Gas hinge for lid drop protection.
- Safety protections are made for the over speed etc to ensure the man machine safety.



WHY CHOOSE US

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific needs.
- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your orders.
- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



Our Factory Introduction

We have a standardized production workshop operated by a dedicated team that ensures quality in every laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production,

quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.
- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.




MKE Exhibition




Packaging & Delivery

MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



After-Sale WARRANTY



FAQ

Q: How long is the production lead time?

A: Within 7-15 workdays after receiving deposit.

Q: What's your warranty?

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

Q: How about the sample of MKE Centrifuge?

A: We could send you the samples but the freight and the sample are charged.

Q: Which kinds of Certificate would you have?

A: ISO9001, ISO13485, CE

Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

Q: Are there customized options available for specific applications?

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



Professional
Centrifuge
Manufacturer

Hunan Michael Laboratory Instrument Co., Ltd.



+86-17775883371



Shirley@mke-lab.com



mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone,
Hunan, China