

ODM Precise Centrifuge Refrigerator Table Top For Veterinary Clinic

Our Product Introduction

for more products please visit us on mkecentrifuge.com

Basic Information

- Certification: CE,ISO13485,ISO9001,SFDA
- Model Number: TRHI22-L (old model:VCK-20R)
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Wooden Box
- Delivery Time: 7-15 work days
- Payment Terms: L/C,D/A,D/P,T/T,Western Union,MoneyGram
- Supply Ability: 100 set/month



Product Specification

- Applications: Hospitals,Veterinary Clinics,Food Labs,Universities & Colleges,Pharmaceutical Labs,Cell Culture Labs
- Max Speed: 22000rpm With Fixed Angle Rotor
12x1.5/2.2ml
- Max RCF: 34090xg
- Max Capacity: 4*250ml Fixed Angle Rotor, 4*750ml Swing Out Rotor
- Temperature Control Range: -20°C +40°C
- Matched Rotors: Fixed Angle Rotors
- Display: LCD
- Warranty: 1 Year Warranty For The Whole Machine,3 Years Warranty,5 Years Warranty,Free Replacement Parts And Shipping Within Warranty
- Highlight: ODM centrifuge refrigerator, Table Top centrifuge refrigerator,



More Images



Product Description

The centrifuge TRHI22-L has an extraordinary capacity with a maximum of 4x250ml fixed angle rotor, and 4x750ml swing out rotor.

A broad selection of fixed-angle and swing-out rotors is available, along with an extensive range of accessories with buckets and adapters, including bio-safety configurations.

The tube capacities extend from 0.2 ml to 250 ml in fixed-angle rotors and 2 ml to 750 ml in swing-out rotors. Microtiter and deepwell plates can also be centrifuged.

The resulting compatibility with a large number of specimen tubes enables high versatility and allows use for a very wide variety of tasks, including high-throughput applications.

Design

Our Product

- Metal housing, metal lid, and stainless steel centrifuge chamber are autoclave-able and resistant detergents and disinfection liquids.
- Observations window on lid for speed calibration and noise measurement.
- Microprocessor control with frequency conversion A.C. motor drive
- 7-inch touch screen display the parameters of Speed, Timer, RCF, Rotors, Program(Memory), Acceleration& Deceleration information and Error code information.
- Double layers sealing ring made for the centrifuge lid and the instrument to keep long good sealing.
- Bio-safety: gas-tight rotating head that can be sterilized at high temperature, centrifugal cavity with ultraviolet sterilization function, enter the sterilization state after turning on the machine.



User friendly

- CFC-free refrigerant, temperature sensor inside the chamber, pre-cooling feature, precise temperature control up to -20°C .
- Programmable and storable, stores up to 1000 programs and provides direct user access with one keystroke save. The stored parameter: combination remains in the memory even after the centrifugation has been switched off.
- Traceability function: USB interface, traceable to more than 1000 sets of centrifugal data.
- Password lock: The centrifuge operating system is equipped with a password to prevent misoperation.
- Provides the 10 acceleration rates and 10 braking ramp, to effectively avoid secondary sedimentation, with better centrifugation effect.
- Automatic lid-locking and holding during rotor run
- Easily exchangeable rotors, special rotor connector which makes it easy to load and unload rotor.
- Automatic troubleshooting with indication of error code and alarm on the screen.
- "Mute"patent technology. Noise level lower than 58 decibel at the highest rotational speed.
- Automatic rpm/rcf conversion
- Automatic rotor identification system, coded rotor identification prevents the rotor from over speeding.



Safety

- Safety and Standards Certification: CE
- Product Safety: EN61010-1 and EN 61010-2-020 Electromagnetic Compatibility: EN 61326-1
- Safety Electrical lid interlock,door can not open when rotor is running and machine won't run when lid open and error code appears.
- Emergency lid-lock release when power failure(very useful in an unexpected power interruption).
- Gas hinge for lid drop protection.
- Safety protections are made for the over speed, over heating(motor, chamber), imbalance etc to ensure the man machine safety.



Hunan Michael Laboratory Instrument Co.,Ltd.

+86-17775883371

Shirley@mke-lab.com

mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone,
Hunan, China