

Portable Microcentrifuge Tabletop Low Speed Centrifuge For PCR Tube

Our Product Introduction

for more products please visit us on mkecentrifuge.com

Basic Information

- Place of Origin: CHINA
- Brand Name: MKE
- Certification: CE, ISO13485, ISO9001, SFDA
- Model Number: ML-6K
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Wooden Box
- Delivery Time: 7-15 work days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 set/month



Product Specification

- Max Speed: 4000/7000rpm
- Max RCF: 2910xg
- Warranty: 1 Year, 3 Years, 5 Years, 1, 1
- Highlight: Portable micro lab equipment,
PCR Tube micro lab equipment,
PCR Tube tabletop low speed centrifuge



More Images



Product Description

Mini centrifuge ML-6K is designed to be compact, portable, and easy to use, making them ideal for small labs or field work. It is an essential tool for many biological and chemical research applications, which is commonly used for DNA extraction, protein purification, and other molecular biology applications.

- Super silent: using special motor, strong power, super silent, stable operation, noise is generally less than 45db, the highest does not exceed 49db.
- Fast instant acceleration: 90% of the maximum speed is reached within 5 seconds.
- Small vibration: well-designed vibration damping device, the movement amplitude is extremely small.
- Easy to replace the rotor: cleverly designed, the rotor can be replaced without any tools.
- High centrifugal force: the relative centrifugal force can reach up to 2910g/5610g.
- Fast stopping: designed with two stopping modes, free stopping in less than 15 seconds and open cover stopping in only 3 seconds.
- High safety factor: there are three touch switches inside, it takes three switches to close at the same time to run, avoiding irregular operation to affect safety.

Our Product Introduction



WHY CHOOSE US

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific needs.
- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your orders.
- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



Our Factory Introduction

We have a standardized production workshop operated by a dedicated team that ensures quality in every laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production, quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.
- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.




MKE Exhibition




Packaging & Delivery

MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



After-Sale WARRANTY



FAQ

Q: How long is the production lead time?

A: Within 7-15 workdays after receiving deposit.

Q: What's your warranty?

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

Q: How about the sample of MKE Centrifuge?

A: We could send you the samples but the freight and the sample are charged.

Q: Which kinds of Certificate would you have?

A: ISO9001, ISO13485, CE

Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

Q: Are there customized options available for specific applications?

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



Professional
Centrifuge
Manufacturer

Hunan Michael Laboratory Instrument Co., Ltd.

+86-17775883371

Shirley@mke-lab.com

mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone, Hunan, China