



Intelligent High Speed Vortex Mixer Shaker Multi Tube For Liquid Powder

Our Product Introduction

for more products please visit us on mkecentrifuge.com

Basic Information

- Place of Origin: China
- Brand Name: MKE
- Certification: CE,ISO13485,ISO9001:2015
- Model Number: VM-8 (old model K8)
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Wooden Box
- Delivery Time: 7-15 work days
- Payment Terms: L/C,D/A,D/P,T/T, Western Union, MoneyGram
- Supply Ability: 100 set/month



Product Specification

- Applications: Liquid, Powder, And Solid Mixing
- Warranty: 1 Year Warranty For The Whole Machine.
- Highlight: Intelligent vortex mixer shaker, Multi Tube vortex mixer shaker, Liquid Powder multi tube vortex mixer



More Images



Product Description

The K8 vortex mixer is a multifunctional mixer designed for quickly mixing liquid, liquid-solid, and powder substances. Its compact and solid structure, as well as its advanced design and manufacturing technology, ensure a simple, smooth and beautiful appearance. The mixer is equipped with the ST intelligent analysis and processing platform, which allows for intelligent control of speed, time, and work mode.

The permanent magnet DC motor is designed with over-temperature protection and no carbon brush, allowing for high-speed mixing without failure or interference. The mixer features two safe operation modes, "point" and "continuous" mode, and a timer function to make experiments more convenient and accurate. The optimized body structure provides a stable and quiet working platform for various oscillatory mixing tasks.

- Minimum Order Quantity: 1
- Price: To be negotiated
- Packaging Details: Plywood Box
- Delivery Time: 7-10 working days
- Payment Terms: T/T
- Supply Ability: 300sets/month

- The K8 vortex mixer is a versatile, easy-to-operate mixer designed to quickly handle liquid, liquid-solid,

Our Product Introduction

and powder mixing.

- Its advanced design and manufacturing ensure a compact, solid structure and a simple, smooth, and beautiful appearance.
- 1.ST intelligent analysis and processing platform can respond to various modes of work requirements, fully equipped with fast processing chipset and electronic sensing technology, easy to achieve intelligent control of speed, time and work mode..
- Permanent magnet DC motor with less heat. Over temperature protection ensures high speed mixing without faults for long time. Brushless design with little spare parts to reduce the operation resistance, noise and interference..
- Two safety operation modes “touch”and “continuous”. Easy switch over. Time count and countdown available for the test convenience.
- Fast thorough mix of any liquid and powder with high efficiency.
- The perfect body structure optimization design, so that it has enough counterweight weight and eccentric balance of operation, for a variety of oscillatory mixing provides a stable, quiet working platform.



Configuration:

K8-1	placing test tube 50mm small container, using slot structure, to prevent the work of the container outside the slip
standard	
K8-2	used to place 85mm flask container class or multiple test tubes mixed together
standard	
K8-3	oscillation diameter 10mm, 14 holes number of test tubes
standard	

Quick Detail:

ST intelligent analysis and processing platform allows for intelligent control of speed, time, and work mode.

Permanent magnet DC motor with over-temperature protection ensures long-term high-speed mixing capabilities with minimal noise.

Two safe operation modes (point and continuous), digital timing, and countdown timer for accurate experimentation.

High speed, high efficiency, and thorough mixing completion in the shortest possible time. Equipped with a stable and quiet working platform for various oscillatory mixing.

WHY CHOOSE US

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific needs.
- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your orders.
- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



Our Factory Introduction

We have a standardized production workshop operated by a dedicated team that ensures quality in every laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production, quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.
- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.




MKE Exhibition




Packaging & Delivery

MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



After-Sale WARRANTY



1 Year
Warranty

Lifetime
maintenance

24-hour
service

15 mins quick
response

Spare parts
lifetime warranty

FAQ

Q: How long is the production lead time?

A: Within 7-15 workdays after receiving deposit.

Q: What's your warranty?

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

Q: How about the sample of MKE Centrifuge?

A: We could send you the samples but the freight and the sample are charged.

Q: Which kinds of Certificate would you have?

A: ISO9001, ISO13485, CE

Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

Q: Are there customized options available for specific applications?

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



Professional
Centrifuge
Manufacturer

Hunan Michael Laboratory Instrument Co.,Ltd.

+86-17775883371

Shirley@mke-lab.com

mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone,
Hunan, China