

ISO9001 Multi Tube Vortex Mixer Lab Equipment For Food Safety Test

Our Product Introduction

for more products please visit us on mkecentrifuge.com


Basic Information

- Place of Origin: China
- Brand Name: MKE
- Certification: CE, ISO13485,ISO9001:2015
- Model Number: MDH12
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Wooden Box
- Delivery Time: 7-15 work days
- Payment Terms: L/C,D/A,D/P,T/T,Western Union,MoneyGram
- Supply Ability: 100 set/month



Product Specification

- Applications: Food Safety Test Of Pesticide Residue, Veterinary Drug Residue, Food Additives
- Warranty: 1 Year Warranty For The Whole Machine.
- Highlight: **ISO9001 multi tube vortex mixer, Food Safety multi tube vortex mixer, ISO9001 vortex lab equipment**



More Images



Product Description

MDH12 multi-tube vortex mixer is mainly for the food safety test of pesticide residue, veterinary drug residue, food additives, and the drug analysis and R&D. Unmatched advantages in handling big quantities of samples. Time and labor saved. Automatic mixing improves the repeatability of experimental methods and the recovery rate of tested samples.

Minimum Order Quantity:1

Price: To be negotiated

Packaging Details: Plywood Box

Delivery Time: 7-10 working days

Payment Terms: T/T

Supply Ability: 300sets/month

- One time 50 samples can be mixed. Standard supply is 10 ml x 50 holes, 15 ml x 50 holes, 20 ml x 50 holes. Non-standard size can be ordered for different tube requirements.
- User-friendly program design. Continuous/time setting/ touch 3 operation modes set for more convenience. Operation ends with beeping sound.
- Perfect structure design ensures sufficient counterweight and circular eccentric balanced operation. Stable and quiet working platform for different samples mixing.



Parts:

Item		Sizes	Note
MDH12-1		10mlx50 hole(diameter 10mm)	Standard
MDH12-2		15mlx50hole (diameter 15mm)	Standard

MDH12-3



20mlx50hole (diameter 20mm) Standard

WHY CHOOSE US

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific needs.
- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your orders.
- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



Our Factory Introduction

We have a standardized production workshop operated by a dedicated team that ensures quality in every laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production, quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.
- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.




MKE Exhibition




Packaging & Delivery

MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



After-Sale WARRANTY



FAQ

Q: How long is the production lead time?

A: Within 7-15 workdays after receiving deposit.

Q: What's your warranty?

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

Q: How about the sample of MKE Centrifuge?

A: We could send you the samples but the freight and the sample are charged.

Q: Which kinds of Certificate would you have?

A: ISO9001, ISO13485, CE

Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

Q: Are there customized options available for specific applications?

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



Professional
Centrifuge
Manufacturer

Hunan Michael Laboratory Instrument Co.,Ltd.



+86-17775883371



Shirley@mke-lab.com



mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone,
Hunan, China