



## Autologous Serum Heating Biofiller Plasma Gel Machine Incubator Thermal Coagulation

Our Product Introduction

for more products please visit us on mkecentrifuge.com


### Basic Information

- Place of Origin: China
- Brand Name: MKE
- Certification: CE, ISO13485, ISO9001:2015
- Model Number: PGM-1 (old model FU-20)
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Wooden Box
- Delivery Time: 7-15 work days
- Payment Terms: L/C,D/A,D/P,T/T,Western Union,MoneyGram
- Supply Ability: 100 set/month



### Product Specification

- Applications: Medical Aesthetics, Mainly For Wrinkle Correction, Body Augmentation Surgery, Filling Of Depressed Scars Caused By Trauma/acne, Etc.
- Warranty: 1 Year Warranty For The Whole Machine.
- Highlight: Autologous biofiller plasma gel machine, Thermal Coagulation biofiller plasma gel machine, Autologous plasma gel filler machine



### More Images



### Product Description

The FU-20 Autologous Serum Heating Incubator (Thermocoagulator) is widely used in the field of medical aesthetics, mainly for wrinkle correction, body augmentation surgery, filling of depressed scars caused by trauma/acne, etc. The Autologous Serum Heating Incubator (Thermocoagulator) is the core equipment for the production of autologous collagen from plasma during the PPP procedure and has a heating function that heats the plasma into a protoplasmic gel state.

Minimum Order Quantity:1  
 Price: To be negotiated  
 Packaging Details: Plywood Box  
 Delivery Time: 5-7 working days  
 Payment Terms: T/T  
 Supply Ability: 500sets/month

- Digital display, temperature and time can be set individually.
- 8 syringes or prp tubes of 10ml-5ml can be heated at the same time, efficient and fast.
- High safety factor, the temperature reaches the set value after the start of the time, time heating is complete sound prompt.

Our Product Introduction



## WHY CHOOSE US

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific needs.
- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your orders.
- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



## Our Factory Introduction

We have a standardized production workshop operated by a dedicated team that ensures quality in every laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production, quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.
- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.




## MKE Exhibition




### **Packaging & Delivery**

MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



### **After-Sale WARRANTY**



### **FAQ**

**Q: How long is the production lead time?**

A: Within 7-15 workdays after receiving deposit.

**Q: What's your warranty?**

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

**Q: How about the sample of MKE Centrifuge?**

A: We could send you the samples but the freight and the sample are charged.

**Q: Which kinds of Certificate would you have?**

A: ISO9001, ISO13485, CE

**Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?**

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

**Q: Are there customized options available for specific applications?**

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



Professional  
Centrifuge  
Manufacturer

Hunan Michael Laboratory Instrument Co., Ltd.

+86-17775883371

Shirley@mke-lab.com

mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone, Hunan, China