


Reagent Card Biochemistry Incubator laboratory equipment For Micro Column Gel Assay Technique

Our Product Introduction

for more products please visit us on mkecentrifuge.com


Basic Information

- Place of Origin: China
- Brand Name: MKE
- Certification: CE, ISO13485, ISO9001:2015
- Model Number: RCI-2(old model FYQ-2)
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Wooden Box
- Delivery Time: 7-15 work days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 set/month




Product Specification

- Applications: Special Instrument Applied To The Micro-column Gel Assay Technique And Suitable For Incubation Of Solution In The Carrier Of Blood Type Reagent Cards At 37°C
- Warranty: 1 Year Warranty For The Whole Machine
- Highlight: Reagent Card biochemistry incubator, Laboratory biochemistry incubator, Reagent Card incubator laboratory equipment



More Images



Product Description


The main purpose of a reagent card incubator is to maintain a consistent and controlled temperature for reagent cards used in diagnostic testing. This is necessary to ensure accurate and reliable results from the testing process.

- Minimum Order Quantity: 1
- Price: To be negotiated
- Packaging Details: Plywood Box
- Delivery Time: 5-7 working days

Our Product

Payment Terms: T/T
Supply Ability: 500sets/month

- LED display, temperature and time parameters can be set.
- With two independent incubation bins, independent timing, complementary impact, each incubation bin can incubate 12 cards at the same time.
- Microprocessor controls temperature and time with good temperature control linearity and small fluctuations.
- It is safe to use, reliable in performance, automatically timed and has an alarm prompt.
- Small size, simple and generous appearance, transparent upper cover, the experimental process is clear at a glance.



WHY CHOOSE US

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific needs.
- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your orders.
- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



Our Factory Introduction

We have a standardized production workshop operated by a dedicated team that ensures quality in every

laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production, quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.
- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.




MKE Exhibition




Packaging & Delivery

MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



After-Sale WARRANTY



FAQ

Q: How long is the production lead time?

A: Within 7-15 workdays after receiving deposit.

Q: What's your warranty?

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

Q: How about the sample of MKE Centrifuge?

A: We could send you the samples but the freight and the sample are charged.

Q: Which kinds of Certificate would you have?

A: ISO9001, ISO13485, CE

Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

Q: Are there customized options available for specific applications?

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



Hunan Michael Laboratory Instrument Co., Ltd.

+86-17775883371

Shirley@mke-lab.com

mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone, Hunan, China