

Rapid Microplate Shaker Incubator For ELISA CLIA Enzyme Linked Immunoassay Incubation

Our Product Introduction

for more products please visit us on mkecentrifuge.com

Basic Information

- Place of Origin: China
- Brand Name: MKE
- Certification: CE,ISO13485,ISO9001:2015
- Model Number: MPI-6(old model FYQ-2)
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Wooden Box
- Delivery Time: 7-15 work days
- Payment Terms: L/C,D/A,D/P,T/T,Western Union,MoneyGram
- Supply Ability: 100 set/month



Product Specification

- Applications: Enzyme-linked Immunosorbent Assay, Incubation Of Various Enzyme Linked Immunoassay (ELISA) And Enzyme Chemiluminescence Immunoassay (CLIA)
- Warranty: 1 Year Warranty For The Whole Machine
- Highlight: **Rapid microplate shaker incubator, CLIA Enzyme microplate shaker incubator, ELISA Enzyme microplate incubator**



Product Description

MK-F6 microplate rapid incubator is a rapid incubator for enzyme-linked immunosorbent assay, which apply special microprocessor temperature control technology. It can be used for the incubation of various enzyme linked immunoassay (ELISA) and enzyme chemiluminescence immunoassay (CLIA).

Minimum Order Quantity:1

Price: To be negotiated

Packaging Details: Plywood Box

Delivery Time: 5-7 working days

Payment Terms: T/T

Supply Ability: 500sets/month

- Rapid incubation in 1/3 the time of conventional incubation. (60 minutes shortened to 20 minutes, 30 minutes shortened to 10 minutes).
- Circulating air heating and humidification, each hole of the enzyme plate is heated evenly.
- Microprocessor control temperature, vibration speed and time can be adjusted, good linear temperature control.
- LCD high definition, large screen display, touch operation interface, clear to understand.
- Vibration speed can be adjusted, soft start running smoothly, to ensure no cross contamination of samples.
- Can Add or take out microplate at any time. Each position can be timed separately according to the set time when operation . There is an audio and video alarm signal at the end of the work, without manual timing.
- Dc brushless motor drive, long life, maintenance free.
- The incubation tray can hold six pcs 96-well microplates at the same time.
- Multiple safety protection functions, in line with the CE safety standard, reminder of fault content and take measures, safe and reliable.



WHY CHOOSE US

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific

needs.

- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your orders.
- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



Our Factory Introduction

We have a standardized production workshop operated by a dedicated team that ensures quality in every laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production, quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.
- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.




MKE Exhibition



Packaging & Delivery

MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



After-Sale WARRANTY



FAQ

Q: How long is the production lead time?

A: Within 7-15 workdays after receiving deposit.

Q: What's your warranty?

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

Q: How about the sample of MKE Centrifuge?

A: We could send you the samples but the freight and the sample are charged.

Q: Which kinds of Certificate would you have?

A: ISO9001, ISO13485, CE

Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

Q: Are there customized options available for specific applications?

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



Professional
Centrifuge
Manufacturer

Hunan Michael Laboratory Instrument Co.,Ltd.



+86-17775883371



Shirley@mke-lab.com



mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone,
Hunan, China