

PPP PRP Biofiller Plasma Gel Maker Machine ISO13485 Certificated

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: MKE
- Certification: CE,ISO13485,ISO9001:2015
- Model Number: tdD-4MC
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Wooden Box
- Delivery Time: 7-15 work days
- Payment Terms: L/C,D/A,D/P,T/T,Western Union,MoneyGram
- Supply Ability: 100 set/month



Product Specification

- Applications: Soften Smile Lines, Restore Lost Fullness In Your Cheeks, Enhance Lips, Contour Your Jaw Line
- Warranty: 1 Year Warranty For The Whole Machine.
- Highlight: **PPP plasma gel maker machine, ISO13485 plasma gel maker machine, PRP biofiller plasma gel machine**



More Images



for more products please visit us on mkecentrifuge.com

Product Description

The TDD-4MC PPP Gel Maker is a quick and reliable device to turn Platelet-Poor Plasma (PPP) into Plasma Bio-Filler Gel.

It can be fabricated by heating at 70 degrees and 90 degrees for 10 minutes. After the production is completed, the machine can automatically cool the prepared PRP gel, which solves the trouble of using the refrigerator to cool. The module is dedicated to 1ml to 5ml syringe.

Minimum Order Quantity:1

Price:	To be negotiated
Packaging Details:	Plywood Box
Delivery Time:	5-7 working days
Payment Terms:	T/T
Supply Ability:	500sets/month

- Two types of machine(digital display & lcd display) for option, the temperature and time can be set separately.
- Equipped with precisely controlled heating and cooling system to turn PPP into Plasma Bio-Filler.
- With two independent chambers: heating chamber and cooling chamber. Separate timing do not affect each other, each chamber can store 10 syringes at the same time, operating fast and efficiently. Safe and reliable performance, automatic timing and alarm prompts.



WHY CHOOSE US

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific needs.
- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your orders.
- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



Our Factory Introduction

We have a standardized production workshop operated by a dedicated team that ensures quality in every laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production, quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.
- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.




MKE Exhibition




Packaging & Delivery

MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



After-Sale WARRANTY



FAQ

Q: How long is the production lead time?

A: Within 7-15 workdays after receiving deposit.

Q: What's your warranty?

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

Q: How about the sample of MKE Centrifuge?

A: We could send you the samples but the freight and the sample are charged.

Q: Which kinds of Certificate would you have?

A: ISO9001, ISO13485, CE

Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

Q: Are there customized options available for specific applications?

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



+86-17775883371



Shirley@mke-lab.com



mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone,
Hunan, China