

## OEM Small Centrifuge Machine Lab Equipment 16500rpm

Our Product Introduction

for more products please visit us on [mkecentrifuge.com](http://mkecentrifuge.com)

### Basic Information

- Place of Origin: CHINA
- Brand Name: MKE
- Certification: CE,ISO13485,ISO9001,SFDA
- Model Number: TG16-WS
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Wooden Box
- Delivery Time: 7-15 work days
- Payment Terms: L/C,D/A,D/P,T/T,Western Union,MoneyGram
- Supply Ability: 100 set/month



### Product Specification

- Applications: Hospitals,Veterinary Clinics,Food Labs,Universities & Colleges,Pharmaceutical Labs,Cell Culture Labs,Infant Stations
- MaxSpeed: 16500rpm
- MaxRCF: 23120xg
- MaxCapacity: 4\*100ml
- Noise Level: ≤58dB
- MatchedRotors: Fixed Angle Rotors
- Display: LCD Display
- Warranty: 1 Year Warranty For The Whole Machine, 3 Years Warranty For The Compressor And 5 Years Warranty For The Motor. Free Replacement Parts And Shipping Within Warranty.
- Highlight: **OEM small centrifuge machine, 16500rpm small centrifuge machine, 16500rpm centrifuge lab equipment**



### More Images



## Product Description

### Description:

Our high-speed centrifuge TG16-WS features a versatile selection of fixed-angle rotors, accommodating PCR tubes, test tubes, and round-bottom tubes. It offers remarkable capacity, holding up to 4 x 100 ml in a fixed angle rotor. With a broad range of rotors and accessories and impressive RCF values reaching up to 23,120 x g, it proves ideal for a diverse array of applications.




### Technical Specifications:

Speed range(rpm)	0-16500rpm adjustable	Speed accuracy	±10rpm
Max Capacity(ml)	4*100ml	RCF range (xg)	0-23120xg adjustable
Display	LCD	Time range	1-99minutes adjustable
Noise(dBA)	≤58dB(A)	Acceleration /Deceleration rates	0-- 10
Dimension(mm)	470x345x295mm	Motor	Frequency Converter motor
Net weight(without rotor)	30KGS	Power supply	AC220V,50HZ,10A ,600W



### Matched Rotors(Customized rotors are available)

No.	Rotor	Capacity	Max.speed	Max.RCF	Note
NO.1 	Angle rotor	12×1.5/2.2ml	16500r/min	19176xg	PP/PC tube
	Adapter	12×0.5/0.2ml			PP/PC tube
NO.2 	Angle rotor	12×5ml	16500r/min	19480xg	PP/PC tube
	Adapter	12×1.5ml			PP/PC tube
NO.3 	Angle rotor	12×10ml	12000r/min	13523xg	PP/PC tube
	Adapter	12×5ml			PP/PC tube
NO.4 	Angle rotor	24×1.5/2.2ml	14000r/min	18407xg	PP/PC tube
	Adapter	24×0.5/0.2ml			PP/PC tube
NO.5 	Angle rotor	36×1.5/2.2ml	13000r/min	15304xg	PP/PC tube
	Adapter	36×0.5/0.2ml			PP/PC tube
NO.6 	Angle rotor	12×8×0.2ml	12000r/min	16985xg	PCR 8-tube strips
NO.7 	Angle rotor	10×15ml	11000r/min	14069xg	PP/PC tube
	Adapter	10×10ml			PP/PC tube
	Adapter	10×5ml			PP/PC tube
NO.8 	Angle rotor	6×50ml	11000r/min	13122xg	PP/PC tube
	Adapter	6×15/10ml			PP/PC tube
	Adapter	6×1.5ml			PP/PC tube
	Angle rotor	4×10/100ml			PP/PC tube

NO.9 	Adapter	4×50ml	10000r/min	10733xg	PP/PC tube
	Adapter	4×15/10ml			PP/PC tube
	Adapter	4×5/1.5ml			PP/PC tube



## Features:

- 1.Compact and space-saving centrifuge
- 2.High Speed Portable Centrifuge TG16-WS featured A.C. Frequency conversion motor direct drive, high precision, low noise, long service life
- 3.Microprocessor control, color break code screen display, program can be programmed and stored, speed, centrifugal force, time can be displayed at the same time.
4. High Speed Portable Centrifuge TG16-WS special Quiet technology, "low noise" patent certificate recognized(patent no. ZL201721835868.9), the highest speed operation noise is less than 58 decibels.
- 5.High Speed Portable Centrifuge 16000rpm TG16-WS features Three levels of damping system, 10 kinds of Lifting and lowering speed curve, 30 self-defined operating procedures can be stored and recalled, can meet different experimental needs.
- 6.With over-speed, over-temperature, electronic safety door lock and other multi-level protection functions to ensure the safety of man and machine.
- 7.Metal housing, metal lid, and stainless steel centrifuge chamber are autoclave-able and resistant detergents and disinfection liquids.

## WHY CHOOSE US

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific needs.
- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your

orders.

- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



### **Our Factory Introduction**

We have a standardized production workshop operated by a dedicated team that ensures quality in every laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production, quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.
- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.




### **MKE Exhibition**






## **Packaging & Delivery**

MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



## **After-Sale WARRANTY**



## **FAQ**

**Q: How long is the production lead time?**

A: Within 7-15 workdays after receiving deposit.

**Q: What's your warranty?**

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

**Q: How about the sample of MKE Centrifuge?**

A: We could send you the samples but the freight and the sample are charged.

**Q: Which kinds of Certificate would you have?**

A: ISO9001, ISO13485, CE

**Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?**

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

**Q: Are there customized options available for specific applications?**

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



Professional  
Centrifuge  
Manufacturer

**Hunan Michael Laboratory Instrument Co.,Ltd.**



+86-17775883371



Shirley@mke-lab.com



mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone,  
Hunan, China