

4000rpm Benchtop Centrifuge Machine Microcentrifuge OEM

Our Product Introduction

Basic Information

- Place of Origin: CHINA
- Brand Name: MKE
- Certification: CE,ISO13485,ISO9001,SFDA
- Model Number: KC-40T
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Wooden Box
- Delivery Time: 7-15 work days
- Payment Terms: L/C,D/A,D/P,T/T,Western Union,MoneyGram
- Supply Ability: 100 set/month




Product Specification

- Applications: Primarily Utilized For Routine Laboratory Tasks, Including The Separation Of Small Volumes Of Biological Samples, Cellular Components, And Various Fluids. These Applications Encompass Processes Like Isolating DNA, RNA, Proteins, And Sedimentation Of Cells
- Max Speed: 4000rpm
- Max RCF: 2680xg
- Max Capacity: 6*50ml Fixed Angle Rotor
- Noise Level: ≤55dB
- Matched Rotors: Fixed Angle Rotors
- Warranty: 1 Year Warranty For The Whole Machine, 3 Years Warranty For The Compressor And 5 Years Warranty For The Motor. Free Replacement Parts And Shipping Within Warranty.
- Highlight: **4000rpm benchtop centrifuge machine , OEM benchtop centrifuge machine**



More Images



for more products please visit us on mkecentrifuge.com

Product Description

Description

Versatility and Precision in One:

Explore unparalleled versatility with our Benchtop Centrifuge KC-40T. Tailored for various applications, it comes with a wide array of fixed-angle rotors, ensuring adaptability to your specific needs.

Accommodating Various Tube Sizes:

Equipped with a diverse selection of accessories, this centrifuge effortlessly accommodates tubes ranging from 2 ml to 50 ml. No matter your sample volume, rest assured, it's handled with precision.

High RCF for Comprehensive Applications:

With a maximum Relative Centrifugal Force (RCF) exceeding 2,680 x g, our KC-40T centrifuge is ideal for a myriad of fields. From clinical diagnostics to veterinary medicine, pharmaceuticals, environmental analysis, and industrial production, it delivers consistent, high-quality results.

Trusted Across Industries and Academia:

Embraced by esteemed universities and research facilities, this centrifuge is a hallmark of reliability and efficiency. Its robust performance ensures it meets the demands of diverse research environments.





Technical Specifications:

Speed range(rpm)	0-4000rpm adjustable	Speed accuracy	±10rpm
Max Capacity(ml)	6*50ml	RCF range (xg)	0-2680xg adjustable
Noise(dBA)	≤55dB(A)	Time range	1-99minutes adjustable
Dimension(mm)	430x320x250mm	Net weight(without rotor)	18KGS

Power supply	AC220V,50HZ,10A, 60W
--------------	----------------------



Matched Rotors

No.	Rotor	Capacity	Max.speed	Max.RCF	Note
NO.1 	Angle rotor	24×10/5ml	4000r/min	2680xg	PP/PC tube / Vacuum blood collection tubes
	Adapter	24×2ml			
	Adapter	24×1.5ml			
NO.2 	Angle rotor	6×50 ml	4000r/min	2680xg	PP/PC tube
	Adapter	6×15/10ml			
	Adapter	6×5ml			



Why Choose Benchtop Centrifuge KC-40T?

Adaptability: Versatile fixed-angle rotors catering to various needs.

Flexible Tube Capacity: Accommodates tubes from 2 ml to 50 ml effortlessly.

Powerful RCF: Achieves an RCF of over 2,680 x g for precise separations.

Wide Application Range: From clinical diagnostics to academia and industrial sectors.

Trusted Performance: Embraced by professionals and researchers globally.

WHY CHOOSE US

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific needs.
- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your orders.
- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



Our Factory Introduction

We have a standardized production workshop operated by a dedicated team that ensures quality in every laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production, quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.
- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.




MKE Exhibition




Packaging & Delivery

MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



After-Sale WARRANTY



FAQ

Q: How long is the production lead time?

A: Within 7-15 workdays after receiving deposit.

Q: What's your warranty?

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

Q: How about the sample of MKE Centrifuge?

A: We could send you the samples but the freight and the sample are charged.

Q: Which kinds of Certificate would you have?

A: ISO9001, ISO13485, CE

Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

Q: Are there customized options available for specific applications?

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



Professional
Centrifuge
Manufacturer

Hunan Michael Laboratory Instrument Co., Ltd.



+86-17775883371



Shirley@mke-lab.com



mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone,
Hunan, China