

## Max.Speed 4000rpm Benchtop Centrifuge Machine Low Speed Separation For Isolating DNA RNA Proteins

Our Product Introduction

for more products please visit us on [mkecentrifuge.com](http://mkecentrifuge.com)

### Basic Information

- Place of Origin: CHINA
- Brand Name: MKE
- Certification: CE,ISO13485,ISO9001,SFDA
- Model Number: TDZ4-WS
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Wooden Box
- Delivery Time: 7-15 work days
- Payment Terms: L/C,D/A,D/P,T/T,Western Union,MoneyGram
- Supply Ability: 100 set/month



### Product Specification

- Applications: Hospitals,Veterinary Clinics,Food Labs,Universities & Colleges,Pharmaceutical Labs,Cell Culture Labs,Infant Stations
- Max Speed: 4000 Rpm
- Max RCF: 2680xg
- Max Capacity: 12\*20ml
- Time Range: 1min-99min Adjustable
- Matched Rotors: Fixed Angle Rotors
- Display: LED
- Noise(dBA): ≤55dB(A)
- Dimension(mm): 430x320x250mm
- Warranty: 1 Year Warranty For The Whole Machine, 3 Years Warranty For The Compressor And 5 Years Warranty For The Motor. Free Replacement Parts And Shipping Within Warranty.
- Highlight: **Benchtop Centrifuge For Isolating ,**



### More Images



## Product Description

### Introducing Benchtop Centrifuge TDZ4-WS:

#### Versatility Redefined:

Meet the Benchtop Centrifuge TDZ4-WS – your versatile lab companion. Equipped with a variety of fixed-angle rotors, it adapts seamlessly to your diverse needs.

#### Accommodating Every Sample:

Its wide range of accessories ensures a perfect fit for tubes spanning from 2 ml to 50 ml. Whether you're handling blood samples in clinical tests, operating in hospital laboratories, or even in beauty salons and dental hospitals, TDZ4-WS delivers consistent results.

#### Precision Where it Matters:

In clinical tests, the accuracy of your results is paramount. TDZ4-WS ensures precise separation, guaranteeing the integrity of your samples and the reliability of your analyses.

#### Across Industries, Trusted Performance:

From hospitals and clinical setups to beauty salons and dental practices, TDZ4-WS is trusted for its reliable performance. Its robust design and efficient operation make it the go-to choice for professionals across diverse fields.



### Technical Specifications:

<b>Speed range(rpm)</b>	0-4000rpm adjustable	<b>Speed accuracy</b>	±10rpm
<b>Max Capacity(ml)</b>	12*20ml	<b>RCF range (xg)</b>	0-2680xg adjustable
<b>Noise(dBA)</b>	≤55dB(A)	<b>Time range</b>	1min-99min adjustable
<b>Dimension(mm)</b>	430x320x250mm	<b>Net weight(without rotor)</b>	18KGS
<b>Power supply</b>	AC220V,50HZ,10A,60W		



## Matched Rotors

No.	Rotor	Capacity	Max.speed	Max.RCF	Note
NO.1	Angle rotor	12×10/5ml	4000r/min	2680xg	PP/PC tube / Vacuum blood collection tubes
	Adapter	12×2ml			Vacuum blood collection tubes
NO.2	Angle rotor	12×20/5ml	4000r/min	2680xg	Vacuum blood collection tubes
	Adapter	12×2ml			Vacuum blood collection tubes
NO.3	Angle rotor	18×10/5ml	4000r/min	2680xg	PP/PC tube/Vacuum blood collection tubes
	Adapter	18×2ml			Vacuum blood collection tubes
	Adapter	18×1.5ml			PP/PC tube
NO.4	Angle rotor	4×50ml	4000r/min	2680xg	PP/PC tube
	Adapter	4×15/10ml			Vacuum blood collection tubes
	Adapter	4×5/2ml			Vacuum blood collection tubes



## Why Choose Benchtop Centrifuge TDZ4-WS?

Adaptable Rotors: Multiple fixed-angle rotors for versatile applications.

Flexible Tube Compatibility: Accommodates tubes from 2 ml to 50 ml seamlessly.

Clinical Precision: Ensures accurate separation for clinical and beauty applications.

Trusted Across Fields: Widely used in hospitals, beauty salons, and dental practices.

Experience the reliability and precision of TDZ4-WS. Where versatility meets accuracy, your lab work finds its perfect partner.

## **WHY CHOOSE US**

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific needs.
- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your orders.
- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



## **Our Factory Introduction**

We have a standardized production workshop operated by a dedicated team that ensures quality in every laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production, quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.
- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.






## **MKE Exhibition**



## **Packaging & Delivery**

MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



## **After-Sale WARRANTY**



## **FAQ**

**Q: How long is the production lead time?**

A: Within 7-15 workdays after receiving deposit.

**Q: What's your warranty?**

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

**Q: How about the sample of MKE Centrifuge?**

A: We could send you the samples but the freight and the sample are charged.

**Q: Which kinds of Certificate would you have?**

A: ISO9001, ISO13485, CE

**Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?**

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

**Q: Are there customized options available for specific applications?**

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



Professional  
Centrifuge  
Manufacturer

**Hunan Michael Laboratory Instrument Co.,Ltd.**



+86-17775883371



Shirley@mke-lab.com



mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone,  
Hunan, China