

4*100ml Benchtop High Speed Laboratory Centrifuge With Ventilated Features And Multiple Rotors

Our Product Introduction

for more products please visit us on mkecentrifuge.com

Basic Information

- Place of Origin: China
- Brand Name: MKE
- Certification: CE, ISO, IUUV
- Model Number: TG16-WS




Product Specification

- Product Type: Benchtop Laboratory Centrifuge
- Matched rotors: Fixed Angle Rotors
- Maximum Capacity: 100 ML X 4
- Max Speed: 16500rpm
- Acceleration /Deceleration Rates: 10 Levels
- Noise Level: ≤58dB
- Timer: 1 Min To 99H 99 Min
- Max Rcf: 19480 Xg
- Display: LED Display
- Highlight: Benchtop High Speed Laboratory Centrifuge, Multiple Rotors Laboratory Centrifuge, Ventilated High Speed Laboratory Centrifuge



More Images



Product Description

4*100ml Benchtop High Speed Laboratory Centrifuge with Ventilated Features and Multiple Rotors **TG16-WS**

	16500 rpm	19480 xg	4x100 ml	≤58 dB(A)
	Maximum Speed	Maximum RCF	Maximum Capacity	Noise Level(dBA)

Technical Specifications:

Speed range(rpm)	0-16500rpm adjustable	Speed accuracy	±10rpm
Max Capacity(ml)	4*100ml	RCF range (xg)	0-19480xg adjustable
Display	LCD	Time range	1-99minutes adjustable
Noise(dBA)	≤58dB(A)	Acceleration /Deceleration rates	0-- 10
Dimension(mm)	470x345x295mm	Motor	Frequency Converter motor
Net weight(without rotor)	30KGS	Power supply	AC220V,50HZ,10A,600W

Matched Rotors(Customized rotors are available)

No.	Rotor	Capacity	Max.speed	Max.RCF	Note
NO.1	Angle rotor	12×1.5/2.2ml	16500r/min	19176xg	PP/PC tube
	Adapter	12×0.5/0.2ml			PP/PC tube
NO.2	Angle rotor	12×5ml	16500r/min	19480xg	PP/PC tube
	Adapter	12×1.5ml			PP/PC tube
NO.3	Angle rotor	12×10ml	12000r/min	13523xg	PP/PC tube
	Adapter	12×5ml			PP/PC tube
NO.4	Angle rotor	24×1.5/2.2ml	14000r/min	18407xg	PP/PC tube
	Adapter	24×0.5/0.2ml			PP/PC tube
NO.5	Angle rotor	36×1.5/2.2ml	13000r/min	15304xg	PP/PC tube
	Adapter	36×0.5/0.2ml			PP/PC tube
NO.6	Angle rotor	12×8×0.2ml	12000r/min	16985xg	PCR 8-tube strips
NO.7	Angle rotor	10×15ml	11000r/min	14069xg	PP/PC tube
	Adapter	10×10ml			PP/PC tube
	Adapter	10×5ml			PP/PC tube
NO.8	Angle rotor	6×50ml	11000r/min	13122xg	PP/PC tube
	Adapter	6×15/10ml			PP/PC tube
	Adapter	6×1.5ml			PP/PC tube
NO.9	Angle rotor	4×10/100ml	10000r/min	10979xg	PP/PC tube
	Adapter	4×50ml			PP/PC tube
	Adapter	4×15/10ml			PP/PC tube
	Adapter	4×5/1.5ml			PP/PC tube

Description:

Our high-speed centrifuge TG16-WS features a versatile selection of fixed-angle rotors, accommodating PCR tubes, test tubes, and round-bottom tubes. It offers remarkable capacity, holding up to 4 x 100 ml in a fixed angle rotor. With a broad range of rotors and accessories and impressive RCF values reaching up to 19480 x g, it proves ideal for a diverse array of applications.

Features

- Compact and space-saving centrifuge
- High Speed Portable Centrifuge TG16-WS featured A.C. Frequency conversion motor direct drive, high precision, low noise, long service life

- Microprocessor control, color break code screen display, program can be programmed and stored, speed, centrifugal force, time can be displayed at the same time.
- High Speed Portable Centrifuge TG16-WS special Quiet technology, "low noise" patent certificate recognized(patent no. ZL201721835868.9), the highest speed operation noise is less than 58 decibels.
- High Speed Portable Centrifuge 16000rpm TG16-WS features Three levels of damping system, 10 kinds of Lifting and lowering speed curve, 30 self-defined operating procedures can be stored and recalled, can meet different experimental needs.
- With over-speed, over-temperature, electronic safety door lock and other multi-level protection functions to ensure the safety of man and machine.
- Metal housing, metal lid, and stainless steel centrifuge chamber are autoclave-able and resistant detergents and disinfection liquids.

WHY CHOOSE US

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific needs.
- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your orders.
- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



Our Factory Introduction

We have a standardized production workshop operated by a dedicated team that ensures quality in every laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production, quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.
- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.




MKE Exhibition




Packaging & Delivery

MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



After-Sale WARRANTY



FAQ

Q: How long is the production lead time?

A: Within 7-15 workdays after receiving deposit.

Q: What's your warranty?

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

Q: How about the sample of MKE Centrifuge?

A: We could send you the samples but the freight and the sample are charged.

Q: Which kinds of Certificate would you have?

A: ISO9001, ISO13485, CE

Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

Q: Are there customized options available for specific applications?

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



Professional
Centrifuge
Manufacturer

Hunan Michael Laboratory Instrument Co.,Ltd.



+86-17775883371



Shirley@mke-lab.com



mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone,
Hunan, China