

4000rpm Speed Adjustable Lab Centrifuge with Max 6 Seats 50ml Capacity

Our Product Introduction

for more products please visit us on mkecentrifuge.com

Basic Information

- Place of Origin: China
- Brand Name: MKE
- Certification: CE, ISO, IUUV
- Model Number: KC-40T




Product Specification

- Type: Benchtop Low Speed Lab Centrifuge
- Timer: 1 Min To 99 Min, Continuous Mode
- Max Capacity: 6*50ml
- Matched Rotors: Fixed Angle Rotors
- Max Speed: 4000rpm
- Safety Features: Lid Lock, Imbalance Detection, Automatic Shut-off
- Power Supply: 110V/220V, 50/60Hz
- Noise Level: ≤58dB
- Highlight: **Max 6 Seats Lab Centrifuge, 50ml Capacity Lab Centrifuge, 4000rpm Lab Centrifuge**




More Images



Product Description

Benchtop Centrifuge KC-40T: Versatility and Precision for Every Lab

Designed for diverse applications, the KC-40T benchtop centrifuge offers unmatched adaptability with its range of fixed-angle rotors and accessories.

	4000 rpm	Maximum Speed
	2268xg	Maximum RCF
	Angle rotor:6x50ml	Maximum Capacity
	≤55dB(A)	Noise Level(dBA)

Accommodates Various Tube Sizes:

From 2 ml to 50 ml, this centrifuge handles all sample volumes with precision and efficiency.

Powerful RCF for Broad Applications:

With a maximum RCF exceeding 2,680 x g, it delivers reliable results across clinical diagnostics, veterinary medicine, pharmaceuticals, environmental studies, and industrial production.

Preferred Across Industries and Academia:

Trusted by top universities and research institutions, the KC-40T ensures robust performance in demanding environments.

Why Choose KC-40T?

Unmatched Versatility: Fixed-angle rotors for diverse applications.

Flexible Capacity: Supports tubes from 2 ml to 50 ml.

Powerful Performance: Achieves RCF over 2,680 x g.

Broad Utility: Ideal for clinical, academic, and industrial use.

Global Trust: Recognized by researchers and professionals worldwide.

Technical Specifications:

Speed range(rpm)	0-4000rpm adjustable	Speed accuracy	±10rpm
Max Capacity(ml)	6*50ml	RCF range (xg)	0-2268xg adjustable
Noise(dBA)	≤55dB(A)	Time range	1-99minutes adjustable
Dimension(mm)	430x320x250mm	Net weight(without rotor)	18KGS
Power supply	AC220V,50HZ,10A, 60W		

Matched Rotors

No.	Rotor	Capacity	Max.speed	Max.RCF	Note
NO.1	Angle rotor	24×10/5ml	4000r/min	2268xg	PP/PC tube / Vacuum blood collection tubes
	Adapter	24×2ml			
	Adapter	24×1.5ml			
NO.2	Angle rotor	6×50 ml	4000r/min	2268xg	PP/PC tube
	Adapter	6×15/10ml			
	Adapter	6×5ml			

WHY CHOOSE US

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific needs.
- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your orders.
- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



Our Factory Introduction

We have a standardized production workshop operated by a dedicated team that ensures quality in every laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production, quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.
- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.



MKE Exhibition



Packaging & Delivery

MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



After-Sale WARRANTY



FAQ

Q: How long is the production lead time?

A: Within 7-15 workdays after receiving deposit.

Q: What's your warranty?

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

Q: How about the sample of MKE Centrifuge?

A: We could send you the samples but the freight and the sample are charged.

Q: Which kinds of Certificate would you have?

A: ISO9001, ISO13485, CE

Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

Q: Are there customized options available for specific applications?

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



Professional
Centrifuge
Manufacturer

Hunan Michael Laboratory Instrument Co., Ltd.



+86-17775883371



Shirley@mke-lab.com



mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone, Hunan, China

