

## Angel Rotor 50ml Tubes Laboratory Centrifuge Small Size RPM RCF Adjustable

Our Product Introduction

for more products please visit us on [mkecentrifuge.com](http://mkecentrifuge.com)

### Basic Information

- Place of Origin: China
- Brand Name: MKE
- Certification: CE, ISO, IUUV
- Model Number: KC-40T



### Product Specification

- Type: 50ml Tubes Laboratory Centrifuge
- Timer: 1 Min To 99 Min, Continuous Mode
- Max Capacity: 6\*50ml
- Matched Rotors: Fixed Angle Rotors
- Max Speed: 4000rpm
- Safety Features: Lid Lock, Imbalance Detection, Automatic Shut-off
- Power Supply: 110V/220V, 50/60Hz
- Noise Level: ≤58dB
- Highlight: 50ml Tubes Laboratory Centrifuge, Small Size Laboratory Centrifuge, Adjustable Laboratory Centrifuge



### More Images



## Product Description

The KC-40T delivers exceptional adaptability with fixed-angle rotors and accessories designed for a variety of lab tasks.

**Supports Multiple Tube Sizes:**

Effortlessly accommodates tubes ranging from 2 ml to 50 ml, ensuring precise performance for any sample size.

**High Centrifugal Force:**

Achieving RCFs over 2,680 x g, it is engineered to excel in clinical diagnostics, research labs, veterinary medicine, and industrial applications.

**Trusted by Leading Institutions:**

Universities, research facilities, and industry professionals rely on the KC-40T for its reliability and efficiency.

**Technical Specifications:**

<b>Speed range(rpm)</b>	0-4000rpm adjustable	<b>Speed accuracy</b>	±10rpm
<b>Max Capacity(ml)</b>	6*50ml	<b>RCF range (xg)</b>	0-2268xg adjustable
<b>Noise(dBA)</b>	≤55dB(A)	<b>Time range</b>	1-99minutes adjustable
<b>Dimension(mm)</b>	430x320x250mm	<b>Net weight(without rotor)</b>	18KGS
<b>Power supply</b>	AC220V, 50HZ, 10A, 60W		



**Matched Rotors**

No.	Rotor	Capacity	Max.speed	Max.RCF	Note
NO.1	Angle rotor	24×10/5ml	4000r/min	2268xg	PP/PC tube / Vacuum blood collection tubes
	Adapter	24×2ml			
	Adapter	24×1.5ml			
NO.2	Angle rotor	6×50 ml	4000r/min	2268xg	PP/PC tube
	Adapter	6×15/10ml			
	Adapter	6×5ml			

**Highlights of the 50ml Tube Lab Centrifuge KC-40T:**

**Versatile Rotors:** Adaptable to a range of applications.

**Wide Capacity Range:** Suitable for tubes between 2 ml and 50 ml.

**Powerful Separation:** RCF exceeding 2,680 x g.

**Comprehensive Applications:** Ideal for academic, clinical, and industrial use.  
**Proven Reliability:** A favorite among global research communities.

## **WHY CHOOSE US**

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific needs.
- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your orders.
- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



## **Our Factory Introduction**

We have a standardized production workshop operated by a dedicated team that ensures quality in every laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production, quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.
- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.



## **MKE Exhibition**






## **Packaging & Delivery**

MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



## **After-Sale WARRANTY**



## **FAQ**

**Q: How long is the production lead time?**

A: Within 7-15 workdays after receiving deposit.

**Q: What's your warranty?**

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

**Q: How about the sample of MKE Centrifuge?**

A: We could send you the samples but the freight and the sample are charged.

**Q: Which kinds of Certificate would you have?**

A: ISO9001, ISO13485, CE

**Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?**

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

**Q: Are there customized options available for specific applications?**

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



Professional  
Centrifuge  
Manufacturer

**Hunan Michael Laboratory Instrument Co.,Ltd.**



+86-17775883371



Shirley@mke-lab.com



mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone,  
Hunan, China