

Versatile TDZ4-WS Laboratory Benchtop Centrifuge Supporting Multiple Rotors Including 4x50ml 12*20ml

Our Product Introduction

for more products please visit us on mkecentrifuge.com

Basic Information

- Place of Origin: CHINA
- Brand Name: MKE
- Certification: CE,ISO13485,ISO9001,SFDA
- Model Number: TDZ4-WS
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Wooden Box
- Delivery Time: 7-15 work days
- Payment Terms: L/C,D/A,D/P,T/T,Western Union,MoneyGram
- Supply Ability: 100 set/month



Product Specification

- Product Type: Versatile TDZ4-WS Laboratory Benchtop Centrifuge
- Applications: Hospitals,Veterinary Clinics,Food Labs,Universities & Colleges,Pharmaceutical Labs,Cell Culture Labs,Infant Stations
- Max Speed: 4000 Rpm
- Max RCF: 2680xg
- Max Capacity: 12*20ml
- Time Range: 1min-99min Adjustable
- Matched Rotors: Fixed Angle Rotors
- Display: LED
- Noise(dBA): ≤55dB(A)
- Dimension(mm): 430x320x250mm
- Warranty: 1 Year Warranty For The Whole Machine, 3 Years Warranty For The Compressor And 5 Years Warranty For The Motor. Free Replacement Parts And Shipping Within



More Images



Product Description

Introducing Benchtop Centrifuge TDZ4-WS:

Unmatched Versatility:

The TDZ4-WS Benchtop Centrifuge redefines adaptability, offering a range of fixed-angle rotors to meet the diverse demands of modern laboratories and clinics.

Sample Handling Made Simple:

With accessories for tubes from 2 ml to 50 ml, this centrifuge is ideal for clinical tests, hospital laboratories, and even specialized applications in beauty and dental practices.

Precision at Its Core:

Engineered for accuracy, the TDZ4-WS ensures reliable sample separation, delivering results you can depend on in critical medical or aesthetic procedures.

Cross-Industry Reliability:

The robust and efficient TDZ4-WS centrifuge is trusted by professionals across industries, combining durability with superior performance for consistent outcomes.



Why Choose Benchtop Centrifuge TDZ4-WS?

Adaptable Rotors: Multiple fixed-angle rotors for versatile applications.

Flexible Tube Compatibility: Accommodates tubes from 2 ml to 50 ml seamlessly.

Clinical Precision: Ensures accurate separation for clinical and beauty applications.


Trusted Across Fields: Widely used in hospitals, beauty salons, and dental practices.

Experience the reliability and precision of TDZ4-WS. Where versatility meets accuracy, your lab work finds its perfect partner.

Technical Specifications:

Speed range(rpm)	0-4000rpm adjustable	Speed accuracy	±10rpm
Max Capacity(ml)	12*20ml	RCF range (xg)	0-2680xg adjustable
Noise(dBA)	≤55dB(A)	Time range	1min-99min adjustable
Dimension(mm)	430x320x250mm	Net weight(without rotor)	18KGS
Power supply	AC220V,50HZ,10A,60W		

Matched Rotors

No.	Rotor	Capacity	Max.speed	Max.RCF	Note
NO.1 	Angle rotor	12×10/5ml	4000r/min	2680xg	PP/PC tube / Vacuum blood collection tubes
	Adapter	12×2ml			Vacuum blood collection tubes
NO.2	Angle rotor	12×20/5ml	4000r/min	2680xg	Vacuum blood collection tubes

NO.3	Adapter	12×2ml	4000r/min	2680xg	Vacuum blood collection tubes
	Angle rotor	18×10/5ml			PP/PC tube/Vacuum blood collection tubes
	Adapter	18×2ml			Vacuum blood collection tubes
	Adapter	18×1.5ml			PP/PC tube
NO.4	Angle rotor	4×50ml	4000r/min	2680xg	PP/PC tube
	Adapter	4×15/10ml			Vacuum blood collection tubes
	Adapter	4×5/2ml			Vacuum blood collection tubes



WHY CHOOSE US

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific needs.
- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your orders.
- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



Our Factory Introduction

We have a standardized production workshop operated by a dedicated team that ensures quality in every laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production, quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.
- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.




MKE Exhibition




Packaging & Delivery

MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



After-Sale WARRANTY



FAQ

Q: How long is the production lead time?

A: Within 7-15 workdays after receiving deposit.

Q: What's your warranty?

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

Q: How about the sample of MKE Centrifuge?

A: We could send you the samples but the freight and the sample are charged.

Q: Which kinds of Certificate would you have?

A: ISO9001, ISO13485, CE

Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

Q: Are there customized options available for specific applications?

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



+86-17775883371



Shirley@mke-lab.com



mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone,
Hunan, China