

High Speed Benthop Centrifuge for Tubes and Bottles 0.2ml 5ml 10ml 20ml 250ml 500ml 750ml

Our Product Introduction

for more products please visit us on mkecentrifuge.com

Basic Information

- Place of Origin: CHINA
- Brand Name: MKE
- Certification: CE,ISO13485,ISO9001,SFDA
- Model Number: VL-200B
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Wooden Box
- Delivery Time: 7-15 work days
- Payment Terms: L/C,D/A,D/P,T/T,Western Union,MoneyGram
- Supply Ability: 100 set/month



Product Specification

- Product Type: Benthop Centrifuge High Speed
- Applications: Hospitals,Veterinary Clinics,Food Labs,Universities & Colleges,Pharmaceutical Labs,Cell Culture Labs,Infant Stations
- Max Speed: 21000rpm
- Max RCF: 31554xg
- Max Capacity: 6*100ml Fixed Angle Rotor, 4*750ml Swing Out Rotor
- Matched Rotors: Fixed Angle Rotors,Swing Out Buckets Rotor,Microplate Rotor
- Display: LCD
- Warranty: 1 Year Warranty For The Whole Machine, 3 Years Warranty For The Compressor And 5 Years Warranty For The Motor. Free Replacement Parts And Shipping Within Warranty.
- Highlight: **High Speed Benthop Centrifuge, Benthop Centrifuge for Tubes,**



More Images



Product Description

Technical Specifications:

| | | | |
|---------------------------|-----------------------|----------------------------------|---------------------------|
| Speed range(rpm) | 0-20000rpm adjustable | Speed accuracy | ±10rpm |
| Max Capacity(ml) | 4*750ml | RCF range (xg) | 0-28174xg |
| Display | LCD | Time range | 1-99minutes adjustable |
| Noise Level(dBA) | ≤58dB(A) | Acceleration /Deceleration rates | 1-- 10 |
| Dimension(mm) | 560x440x360mm | Motor | Frequency Converter motor |
| Net weight(without rotor) | 45KGS | Power supply | AC220V,50HZ,10A,770W |

Elevate your lab's capabilities with the VL-200B high-speed centrifuge, designed for maximum efficiency and precision. Ideal for diverse applications, its compact size and quiet performance make it a standout choice.

Matched Rotors(Customized rotors are available)

| No. | Rotor | Capacity | Max. speed | Max.RCF | Note |
|------|-------------|--------------|------------|---------|------------|
| NO.1 | Angle rotor | 12×1.5-2.2ml | 20000r/min | 28174xg | PP/PC tube |
| | Adapter | 12×0.5ml | | | PP/PC tube |
| | Adapter | 12×0.2ml | | | PP/PC tube |
| NO.2 | Angle rotor | 12×5ml | 17000r/min | 20680xg | PP/PC tube |
| | Adapter | 12×1.5ml | | | PP/PC tube |
| NO.3 | Angle rotor | 24×1.5-2.2ml | 16000r/min | 24041xg | PP/PC tube |
| | Adapter | 24×0.5ml | | | PP/PC tube |
| | Adapter | 24×0.2ml | | | PP/PC tube |
| NO.4 | Angle rotor | 10×15ml | 15000r/min | 26161xg | PP/PC tube |
| | Adapter | 10×10ml | | | PP/PC tube |
| | Adapter | 10×5ml | | | PP/PC tube |
| | Adapter | 10×1.5ml | | | PP/PC tube |
| NO.5 | Angle rotor | 48×1.5-2.2ml | 15000r/min | 23394xg | PP/PC tube |
| | Adapter | 48×0.5ml | | | PP/PC tube |
| | Adapter | 48×0.2ml | | | |
| NO.6 | Angle rotor | 12×10ml | 15000r/min | 21130xg | PP/PC tube |
| | Adapter | 12×5ml/1.5ml | | | PP/PC tube |
| NO.7 | Angle rotor | 6×50ml | 14000r/min | 21255xg | PP/PC tube |
| | Adapter | 6×15ml | | | PP/PC tube |
| | Adapter | 6×10ml | | | PP/PC tube |
| | Adapter | 6×1.5ml | | | PP/PC tube |
| NO.8 | Angle rotor | 6×100ml | 11000r/min | 14610xg | PP/PC tube |
| | Adapter | 6×50ml | | | PP/PC tube |
| | Adapter | 6×15ml | | | PP/PC tube |
| | Adapter | 6×10ml | | | PP/PC tube |
| | Adapter | 6×5ml | | | PP/PC tube |
| | Adapter | 6×1.5ml | | | PP/PC tube |

| | | | | | |
|-------|--------------|--------------|------------|---------|------------------------------|
| NO.9 | Angle rotor | 8×50ml | 12000r/min | 17226xg | PP/PC tube |
| | Adapter | 8×15ml | | | PP/PC tube |
| | Adapter | 8×10ml | | | PP/PC tube |
| | Adapter | 8×1.5ml | | | PP/PC tube |
| NO.10 | Swing bucket | 4×250ml | 4000r/min | 3049xg | Round cup |
| | Adapter | 4×9×2-5ml | | | Vacuum blood collection tube |
| | Adapter | 4×5×10-15ml | | | PP/PC tube |
| | Adapter | 4×50ml | | | PP/PC tube |
| | Adapter | 4×100ml | | | PP/PC tube |
| | Adapter | 4×9×1.5ml | | | PP/PC tube |
| NO.11 | Swing bucket | 4×400ml | 4000r/min | 3049xg | Square cup |
| | Adapter | 4×16×2-5ml | | | Vacuum blood collection tube |
| | Adapter | 4×9×10-15ml | | | PP/PC tube |
| | Adapter | 4×4×50ml | | | PP/PC tube |
| | Adapter | 4×100ml | | | PP/PC tube |
| | Adapter | 4×15×1.5ml | | | PP/PC tube |
| NO.12 | Swing bucket | 4×500ml | 4000r/min | 3354xg | Round cup |
| | Adapter | 4×17×2-5ml | | | Vacuum blood collection tube |
| | Adapter | 4×8×10-15ml | | | PP/PC tube |
| | Adapter | 4×3×50ml | | | PP/PC tube |
| | Adapter | 4×100ml | | | PP/PC tube |
| | | | | | |
| NO.13 | Swing bucket | 4×600ml | 4000r/min | 3354xg | Square cup |
| | Adapter | 4×25×2-5ml | | | Vacuum blood collection tube |
| | Adapter | 4×12×10-15ml | | | PP/PC tube |
| | Adapter | 4×5×50ml | | | PP/PC tube |
| | Adapter | 4×4×100ml | | | PP/PC tube |
| | | | | | |
| NO.14 | Swing bucket | 4×750ml | 4000r/min | 3443xg | Round cup |
| | Adapter | 4×26×2-5ml | | | Vacuum blood collection tube |
| | Adapter | 4×13×10-15ml | | | PP/PC tube |
| | Adapter | 4×5×50ml | | | PP/PC tube |

| | | | | | |
|-------|------------------|-------------|-----------|--------|------------|
| | Adapter | 4×3×100ml | | | PP/PC tube |
| NO.15 | Microplate rotor | 2×2×96 well | 4200r/min | 2860xg | Square cup |
| | Adapter | 2×32×1.5ml | | | PP/PC tube |

Key Features:

High Performance: Speeds up to 21,000 rpm, RCF up to 31,554 x g, and tube capacities from 0.2 ml to 750 ml.

Smart Design: Low noise with patented “Mute” technology, automatic rotor recognition, and adjustable acceleration/braking rates.

Easy to Use: Intuitive LCD interface, quick rotor changes, and memory retention for stored settings.

Safe and Reliable: IEC 61010 compliance, emergency lid release, and anti-vibration damping ensure smooth and secure operation.

The VL-200B centrifuge combines advanced functionality with user-friendly features, delivering reliable results every time.

WHY CHOOSE US

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific needs.
- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your orders.
- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



Our Factory Introduction

We have a standardized production workshop operated by a dedicated team that ensures quality in every laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production, quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.

- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.




MKE Exhibition



Packaging & Delivery

MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



After-Sale WARRANTY



FAQ

Q: How long is the production lead time?

A: Within 7-15 workdays after receiving deposit.

Q: What's your warranty?

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

Q: How about the sample of MKE Centrifuge?

A: We could send you the samples but the freight and the sample are charged.

Q: Which kinds of Certificate would you have?

A: ISO9001, ISO13485, CE

Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

Q: Are there customized options available for specific applications?

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



Professional
Centrifuge
Manufacturer

Hunan Michael Laboratory Instrument Co.,Ltd.



+86-17775883371



Shirley@mke-lab.com



mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone,
Hunan, China