

## 2000rpm To 4000 Rpm Low-Speed Centrifuge 15ml 10ml 5ml Test Tubes Laboratory Use

Our Product Introduction

for more products please visit us on [mkecentrifuge.com](http://mkecentrifuge.com)

### Basic Information

- Place of Origin: CHINA
- Brand Name: MKE
- Certification: CE,ISO13485,ISO9001,SFDA
- Model Number: TD4
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Wooden Box
- Delivery Time: 7-15 work days
- Payment Terms: L/C,D/A,D/P,T/T,Western Union,MoneyGram
- Supply Ability: 100 set/month



### Product Specification

- Product Type: 2000rpm To 4000 Rpm Low-Speed Centrifuge
- Applications: Hospitals,Veterinary Clinics,Food Labs,Universities & Colleges,Pharmaceutical Labs,Cell Culture Labs
- Max.Speed: 4000rpm
- Max.RCF: 2200xg
- Max.Capacity: 8\*15ml Angle Rotor
- Matched Rotors: Fixed Angle Rotors
- Display: Digital
- Warranty: 1 Year,Free Parts Within Warranty
- Highlight: 5ml Low Speed Centrifuge,  
10ml Low Speed Centrifuge,  
15ml Low Speed Centrifuge



### More Images



## Product Description

Designed for versatility and ease of use, the TD4 low-speed centrifuge is ideal for hospitals, research labs, and beauty clinics. Its compact size and lightweight construction ensure effortless operation while delivering reliable results.

### Why Choose TD4?

#### Smart and Efficient

Automatic lid-lock and quick rotor exchange simplify tasks.

Real-time r/min to RCF conversion guarantees precision.

Intelligent troubleshooting with on-screen error codes and alerts.

Operates quietly with patented "Mute" technology, maintaining noise levels below 55 dB.

#### Built to Last

Sturdy ABS plastic exterior withstands frequent cleaning with disinfectants.

Aluminum alloy rotor ensures dependable performance for years.

#### Safety First

Compliant with IEC 61010 standards for secure operation.

Features include an electrical lid interlock, emergency lid release, and overspeed protection.

Gas hinges and other safeguards prioritize user safety.

#### User-Centric Design

Digital display for clear monitoring of key parameters.

Observation window allows easy speed and noise assessments.

Compact and lightweight design offers excellent portability.

#### Technical Specifications:

<b>Speed range(rpm)</b>	0-4000rpm adjustable	<b>Speed accuracy</b>	±10rpm
<b>Max Capacity(ml)</b>	8*15ml	<b>RCF range (xg)</b>	0-2200xg adjustable
<b>Noise(dBA)</b>	≤55dB(A)	<b>Time range</b>	0-99minutes adjustable
<b>Dimension(mm)</b>	320x270x215mm	<b>Net weight(without rotor)</b>	7KGS
<b>Power supply</b>	AC220V,50HZ,10A,60W		

#### Matched Rotors

No.	Rotor	Capacity	Max. speed	Max. RCF	Note
NO.1	Angle rotor	8×15/10/5ml	4000r/min	1900xg	Blood collection tube/PP/PC tube
	Adapter	8×2ml			Blood collection tube
NO.2	Angle Rotor	12×5ml	4000r/min	2100xg	Blood collection tube
	Adapter	12×2ml			Blood collection tube

### WHY CHOOSE US

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific needs.
- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your orders.
- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



## Our Factory Introduction

We have a standardized production workshop operated by a dedicated team that ensures quality in every laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production, quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.
- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.




## MKE Exhibition




## Packaging & Delivery



MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



### **After-Sale WARRANTY**



### **FAQ**

**Q: How long is the production lead time?**

A: Within 7-15 workdays after receiving deposit.

**Q: What's your warranty?**

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

**Q: How about the sample of MKE Centrifuge?**

A: We could send you the samples but the freight and the sample are charged.

**Q: Which kinds of Certificate would you have?**

A: ISO9001, ISO13485, CE

**Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?**

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

**Q: Are there customized options available for specific applications?**

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



+86-17775883371



Shirley@mke-lab.com



mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone,  
Hunan, China