

Low-Speed Laboratory Centrifuge Compact Power For Precision Tasks OEM/ODM Support

Our Product Introduction

for more products please visit us on mkecentrifuge.com

Basic Information

- Place of Origin: CHINA
- Brand Name: MKE
- Certification: CE,ISO13485,ISO9001,SFDA
- Model Number: TD4
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Wooden Box
- Delivery Time: 7-15 work days
- Payment Terms: L/C,D/A,D/P,T/T,Western Union,MoneyGram
- Supply Ability: 100 set/month



Product Specification

- Product Type: Low-Speed Laboratory Centrifuge Compact Power For Precision Tasks
- Applications: Hospitals,Veterinary Clinics,Food Labs,Universities & Colleges,Pharmaceutical Labs,Cell Culture Labs
- Max.Speed: 4000rpm
- Max.RCF: 2200xg
- Max.Capacity: 8*15ml Angle Rotor
- Matched Rotors: Fixed Angle Rotors
- Display: Digital
- Warranty: 1 Year,Free Parts Within Warranty
- Highlight: **Compact Power Laboratory Centrifuge, OEM/ODM Laboratory Centrifuge, Precision Tasks Laboratory Centrifuge**



More Images



Product Description

Technical Specifications:

Speed range(rpm)	0-4000rpm adjustable	Speed accuracy	±10rpm
Max Capacity(ml)	8*15ml	RCF range (xg)	0-2200xg adjustable
Noise(dBA)	≤55dB(A)	Time range	0-99minutes adjustable
Dimension(mm)	320x270x215mm	Net weight(without rotor)	7KGS
Power supply	AC220V,50HZ,10A,60W		

The TD4 centrifuge is the ideal choice for small-volume separations, perfectly suited for hospitals, research labs, and beauty centers. Its lightweight and compact design make it portable and easy to use.

Matched Rotors

No.	Rotor	Capacity	Max. speed	Max. RCF	Note
NO.1	Angle rotor	8×15/10/5ml	4000r/min	1900xg	Blood collection tube/PP/PC tube
	Adapter	8×2ml			Blood collection tube
NO.2	Angle Rotor	12×5ml	4000r/min	2100xg	Blood collection tube
	Adapter	12×2ml			Blood collection tube

Key Features:

Smart Technology: Automatic lid-lock, real-time r/min to RCF conversion, and troubleshooting with error codes.

Quiet Operation: Patented "Mute" technology keeps noise levels below 55 dB.

Robust Design: Durable ABS plastic body and aluminum alloy rotor ensure long-term performance.

Safety Guaranteed: IEC 61010 compliance with features like lid interlocks, emergency lid release, and overspeed protection.

User-Friendly: Digital display, observation window, and compact structure for convenience.

Choose the TD4 centrifuge for reliable, efficient, and safe separation tasks across various fields.

WHY CHOOSE US

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific needs.
- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your orders.
- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



Our Factory Introduction

We have a standardized production workshop operated by a dedicated team that ensures quality in every laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production, quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.
- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.



MKE Exhibition




Packaging & Delivery

MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



After-Sale WARRANTY



FAQ

Q: How long is the production lead time?

A: Within 7-15 workdays after receiving deposit.

Q: What's your warranty?

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

Q: How about the sample of MKE Centrifuge?

A: We could send you the samples but the freight and the sample are charged.

Q: Which kinds of Certificate would you have?

A: ISO9001, ISO13485, CE

Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

Q: Are there customized options available for specific applications?

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



+86-17775883371



Shirley@mke-lab.com



mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone,
Hunan, China