

## Brushless AC Motor Floor-Standing High-Speed Refrigerated Centrifuge 4000 Rpm To 20000 Rpm

Our Product Introduction

for more products please visit us on [mkecentrifuge.com](http://mkecentrifuge.com)


### Basic Information

- Place of Origin: China
- Brand Name: MKE
- Certification: CE, ISO, IUUV
- Model Number: VCK-21R
- Minimum Order Quantity: 1
- Price: Negotiate
- Packaging Details: Plywood Box
- Delivery Time: 10-15 work days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 50 units/month



### Product Specification

- Product: Brushless AC Motor Floor-Standing High-Speed Refrigerated Centrifuge
- Temperature Control Range: -20°C +40°C
- Maximum Capacity: 4x750ml
- Rcf: 26832 Xg
- Acceleration Profiles: 1~10
- Rotor Identification: Automatic
- Certifications: CE, ISO9001
- Timer: 1-99 Minutes
- Max Speed: 20000rpm
- Highlight: Brushless AC Motor Refrigerated Centrifuge, Floor Standing Refrigerated Centrifuge, High Speed Refrigerated Centrifuge



### More Images



## Product Description

The VCK-21R centrifuge is ideal for laboratory applications, combining high-speed performance with versatile functionality and robust safety measures.

### Why Choose the VCK-21R?

#### Efficiency and Power:

High-power AC motor provides precision, reliability, and a long service life.

Advanced refrigeration system ensures accurate temperature control and eco-friendly operation.

#### User-Friendly Design:

Full-color screen displays all critical parameters simultaneously for ease of use.

Supports both high-speed (20,000 rpm) and large-capacity separations (4 x 750 ml).

Noise-reduction technology ensures quiet operation (<58 dB).

#### Customizable and Safe:

Programmable settings with up to 30 stored procedures adapt to different experimental needs.

Adapters accommodate a wide range of tube sizes, from 0.5 ml to 500 ml.

Multi-level safety features protect both users and equipment.

### Technical Specifications:

Speed range(rpm)	0-20000rpm adjustable	Time range	1-99minutes 59s adjustable
Max Capacity(ml)	4*750ml	Speed accuracy	±10r/min
RCF range (xg)	0-26832xg adjustable	ACC/DEC	1-10 rates
Temperature control range	-20°C +40°C	Temperature control accuracy	±1°C
Noise(dBA)	≤58dB(A)	Net weight(without rotor)	190KGS
Power supply	AC220V,50HZ,16A, 2150W	Dimension(mm)	700x600x820mm(LxW xH)

### Matched Rotors

No.	Rotor	Max. capacity	Max.speed	Max.RCF	Note
NO.1	Angle rotor	10×5ml	20000r/min	26832xg	PP/PC tube
	Adapter	10×1.5ml			PP/PC tube
NO.2	Angle rotor	24×1.5/2.2ml	16000r/min	24041xg	PP/PC tube
	Adapter	24×0.5ml			PP/PC tube
	Adapter	24×0.2ml			PP/PC tube
NO.3	Angle rotor	12×10ml	15000r/min	21130xg	PP/PC tube
	Adapter	12×5ml			PP/PC tube
	Adapter	12×1.5ml			PP/PC tube
NO.4	Angle rotor	10×15ml	15000r/min	26161xg	PP/PC tube
	Adapter	10×5ml			PP/PC tube
	Adapter	10×1.5ml			PP/PC tube
NO.5	Angle rotor	6×50ml	15000r/min	24400xg	PP/PC tube
	Adapter	6×15ml			PP/PC tube
	Adapter	6×10ml			PP/PC tube
	Adapter	6×1.5ml			PP/PC tube
NO.6	Angle rotor	8×50ml	12000r/min	17226xg	PP/PC tube
	Adapter	8×15ml			PP/PC tube
	Adapter	8×10ml			PP/PC tube
	Adapter	8×1.5ml			PP/PC tube
NO.7	Angle rotor	6×100ml	13000r/min	20406xg	PP/PC tube
	Adapter	6×50ml			PP/PC tube
	Adapter	6×15ml			PP/PC tube
	Adapter	6×10ml			PP/PC tube
	Adapter	6×5ml			PP/PC tube
	Adapter	6×1.5ml			PP/PC tube
NO.8	Angle rotor	4×250ml	10000r/min	14758xg	PP/PC tube
	Adapter	4×100ml			PP/PC tube
	Adapter	4×50ml			PP/PC tube
	Adapter	4×5×15ml			PP/PC tube
	Adapter	4×5×10ml			PP/PC tube
	Adapter	4×10×1.5ml			PP/PC tube
NO.9	Swing bucket	4×600ml	4000r/min	3354xg	Square cup
	Adapter	4×25×5/2ml			PP/PC tube
	Adapter	4×12×15/10ml			PP/PC tube
	Adapter	4×5×50ml			PP/PC tube
	Adapter	4×4×100ml			PP/PC tube

NO.10	Swing rotor	4×500ml	4000r/min	3354xg	Round cup
	Adapter	4×17×5/2ml			PP/PC tube
	Adapter	4×8×15/10ml			PP/PC tube
	Adapter	4×3×50ml			PP/PC tube
	Adapter	4×100ml			PP/PC tube
NO.11	Swing rotor	4×750ml	4000r/min	3443xg	Round cup
	Adapter	4×26×5/2ml			PP/PC tube
	Adapter	4×13×15/10ml			PP/PC tube
	Adapter	4×5×50ml			PP/PC tube
	Adapter	4×3×100ml			PP/PC tube
NO.12	Microplate rotor	4×2×96 well	4000r/min	2978xg	4*2*96 deep well plate 4*4*96 shallow well plate
	Adapter	4×32×1.5ml			PP/PC tube

## WHY CHOOSE US

- **Direct from the Manufacturer:** Enjoy competitive pricing on all our products.
- **Certified Excellence:** Experience premium quality with our CE and ISO certified products.
- **Customization at Its Best:** Benefit from our robust OEM and ODM capabilities tailored to your specific needs.
- **Reliable Production and Delivery:** Trust us for guaranteed on-time production and delivery of your orders.
- **Comprehensive Support:** Rely on our comprehensive after-sales service, including complimentary product accessories.
- **Expertise You Can Trust:** We stand behind our 15 years of experience in R&D and production of medical centrifuge machines.



## Our Factory Introduction

We have a standardized production workshop operated by a dedicated team that ensures quality in every laboratory centrifuge and lab equipment we produce. Our R&D division includes design engineers, technology specialists, qualified experts, and process control managers. Following a strict quality control system, our specialists oversee the entire production process, from design and R&D to testing, production, quality inspection, and packaging.

- Specialize in centrifuges.
- With decades of experience, we've exported to 80+ countries.
- Proven Reputation: SGS and Intertek inspected, underlining our reliability.
- Manufacturing Capabilities: 6 floors, 3000 sqm each, dedicated to manufacturing and quality control.






## **MKE Exhibition**




## **Packaging & Delivery**

MKE centrifuges are securely packaged in wooden cases or carton boxes, selected according to the weight and size of the goods. This ensures they are delivered in excellent condition to customers worldwide.



## **After-Sale WARRANTY**



## **FAQ**

**Q: How long is the production lead time?**

A: Within 7-15 workdays after receiving deposit.

**Q: What's your warranty?**

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

**Q: How about the sample of MKE Centrifuge?**

A: We could send you the samples but the freight and the sample are charged.

**Q: Which kinds of Certificate would you have?**

A: ISO9001, ISO13485, CE

**Q: Do you offer centrifuges with voltage and plug options suitable for my country/region?**

A: Yes, we provide centrifuges with various voltage and plug options to meet different country/region requirements. Please let us know your specific needs, and we will ensure compatibility.

**Q: Are there customized options available for specific applications?**

A: Absolutely! We offer customization options to meet specific application needs. Let us know your requirements, and we will tailor the centrifuge accordingly.



**Hunan Michael Laboratory Instrument Co.,Ltd.**



+86-17775883371



Shirley@mke-lab.com



mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone,  
Hunan, China