

## Low Speed Benchtop Lab Centrifuge with 4\*600ml Swing Rotor and 6080xg RCF for Hospitals and Universities

Our Product Introduction

for more products please visit us on [mkecentrifuge.com](http://mkecentrifuge.com)

### Basic Information

- Place of Origin: Hunan, China
- Brand Name: MKE
- Certification: CE, ISO13485, ISO9001, SFDA
- Model Number: TA6-L (old model: KC-500B)
- Minimum Order Quantity: 1
- Price: \$892-935



### Product Specification

- Max Speed: 6000 Rpm
- Max RCF: 6080xg
- Max Capacity: 4\*600ml Swing Out Rotor
- Matched Rotors: Fixed Angle Rotors + Swing Out Rotors
- Display Type: LED
- Speed Accuracy:  $\pm 10$ rpm
- Time Range: 1-99h59min59s
- Net Weight: 50KG
- Noise Level:  $\leq 55$ dB(A)
- Acceleration/Deceleration: 0-9
- Dimensions: 500x400x300mm
- Motor Type: Frequency Converter Motor
- Power Supply: AC220V, 50HZ, 5A, 250W
- Material: Metal, Plastic, Steel
- Instrument Classification: Class I



### More Images



## Product Description

### Low Speed Benchtop Lab Centrifuge Machine 4×600ml Swing Rotor Centrifuge with Automatic Off-cap Tec

#### Premium Laboratory Centrifuge with Comprehensive Warranty

**All MKE Centrifuge: One-year full unit warranty, five-year motor warranty**

The T5-L High-Quality Unrefrigerated Low Speed Benchtop Centrifuge combines advanced technology with high capacity. Featuring multiple swing-out and fixed-angle rotors, it accommodates tubes from 2ml to 100ml, including 15ml, 50ml, and vacuum blood collection tubes. With a maximum RCF of  $4,053 \times g$ , it's ideal for diverse applications in clinical diagnostics, veterinary medicine, pharmaceuticals, environmental analysis, and industrial production.

#### Key Specifications

Applications	Hospitals, Veterinary clinics, Food labs, Universities & colleges, Pharmaceutical labs, Cell culture labs, Infant stations
Max. Speed	5000 rpm
Max. RCF	4053×g
Max. Capacity	4×100ml swing out rotor
Matched Rotors	Fixed angle rotors + Swing out rotor
Display	LED
Model Number	T5-L (old model: KC-50X)
Certification	CE, ISO13485, ISO9001, SFDA

#### Warranty Information

1 year warranty for the whole machine, 3 years warranty for the compressor, and 5 years warranty for the motor. Includes free replacement parts and shipping within warranty period.

#### Technical Specifications

Speed range (rpm)	0-5000rpm adjustable	Speed accuracy	±10rpm
Max Capacity (ml)	4×100ml	RCF range (×g)	0-4053×g adjustable
Display	LED	Time range	0-99 minutes 59s
Noise (dBA)	≤55dB(A)	Acceleration/Deceleration rates	0-9
Dimension (mm)	500×400×300mm	Motor	Frequency Converter motor
Net weight (without rotor)	25KGS	Power supply	AC220V, 50HZ, 5A, 250W

#### Available Accessories (All rotors are autoclavable and customizable)

No.	Item	Capacity	Max. speed	Max. RCF	Note
NO.1	Angle rotor	12×15/10ml	5000r/min	3047×g	PP/PC tube
	Adapter	12×5/2ml			Vacuum blood collection tube

Custom adapters are available for specific requirements. Contact us directly for more details about customization options.

#### Our Advantages

### Warranty & Service

1 year for the entire machine, 3 years for the compressor, and 5 years for the motor

Extended warranty option available

Factory maintenance includes freight and customs

Lifelong technical support with remote assistance and on-site training

Distributor assistance for installation and maintenance

### Frequently Asked Questions

**Q: How long is the production lead time?**

A: Within 7-15 workdays after receiving deposit.

**Q: What's your warranty?**

A: All MKE products have ONE year warranty, five-year motor warranty. We won't charge clients any money for repairing the machine during one-year warranty time, the freight is also on us.

**Ready to Order or Need More Information?**

**Contact Now and Get Extra  
Discount**



Professional  
Centrifuge  
Manufacturer

**Hunan Michael Laboratory Instrument Co.,Ltd.**



+86-17775883371



Shirley@mke-lab.com



mkecentrifuge.com

MKE Building, Sanhe Zhihui Industrial Park, Guyuan Road No. 426, Changsha High-tech Development Zone,  
Hunan, China